





# HOW TO USE OPERA'S PUMPING MODES

Vacuum strength and cycle speeds can be adjusted manually and uniquely for either breast, enabling a personalized experience.

 Toggle between EXPRESSION and MASSAGE by pressing this button.

 Press this button to quickly change modes. Single Left, Single Right, Alternate, and Synchronous current pumping mode will be displayed on the top-right corner

## ADJUSTING BOTH SIDES AT THE SAME TIME

When Opera is in **SYNCHRONOUS** or **ALTERNATE** mode, adjusting the vacuum or cycle will increase or decrease both sides at the same time.  
**SYNCHRONOUS** mode express both breasts at the same time.  
**ALTERNATE** mode express both breasts in an alternate pattern.



Increasing or decreasing the vacuum level and/or cycle will affect both sides.



Screen will show vacuum level and cycle for both sides.

## ADJUSTING ONE SIDE

Opera's dual motors technology allows moms to adjust each side independently.  
Press **MODE** to select the side you want to use. (Left or Right)  
Once the side is selected, you can modify the vacuum level and cycle for that side.



Increasing or decreasing the vacuum level and/or cycle will affect the currently selected side only.




Screen will show the current vacuum level and cycle for the selected side.

## USING DIFFERENT SUCTION LEVELS ON EACH SIDE WHILE DOUBLE PUMPING


Once you have adjusted each side to your preference, you can press **MODE** to select **SYNCHRONOUS** or **ALTERNATE** mode.  
Opera will keep the vacuum level selected for each side. However, cycles will be the same for both sides.  
Opera will use the lowest cycle setting for both sides.

### SYNCHRONOUS OR ALTERNATE MODE



Same Cycle  
Different Level

Increasing or decreasing the vacuum level and/or cycle will affect both sides. But the difference in suction level between sides will remain



Once one side reached max vacuum level, if you keep increasing the level it will only affect the side that had a lower level until both sides reach the maximum level.

The same will happen if you decrease the vacuum level. The difference in vacuum will remain until both reach the lowest level. (0)